

**MHC Class II I-Ab Mouse anti-Mouse Monoclonal (FITC) Antibody - LS-C72619 - LSBio**

<b>CatalogID:</b>	LS-C72619
<b>Target:</b>	MHC Class II I-Ab
<b>Host</b>	Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgM,k
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	Mouse
<b>Antigen Type:</b>	Cells
<b>Immunogen:</b>	C3H. SW splenocytes.
<b>Specificity:</b>	Reacts specifically with the I-Ab haplotype of MHC Class II molecules.
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, sodium azide.
<b>Usage Summary:</b>	Flow cytometry. Complement-dependent cytotoxicity assays. Flow Cytometry: 1 ug/1E6 cells.
<b>Uses:</b>	Flow Cytometry, Cytotoxic Assay (Optimal dilution to be determined by the researcher)
<b>Size:</b>	500 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences