

MHC Class II I-Ab Mouse anti-Mouse Monoclonal (PE) Antibody - LS-C72621 - LSBio

CatalogID:	LS-C72621
Target:	MHC Class II I-Ab
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgM,k
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Mouse
Antigen Type:	Cells
Immunogen:	C3H. SW splenocytes.
Specificity:	Reacts specifically with the I-Ab haplotype of MHC Class II molecules.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Usage Summary:	Flow cytometry. Complement-dependent cytotoxicity assays. Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	Flow Cytometry, Cytotoxic Assay (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences