

MHC Class II B-L Mouse anti-Chicken Monoclonal (FITC) Antibody - LS-C72632 - LSBio

CatalogID:	LS-C72632
Target:	MHC Class II B-L
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Chicken
Specificity:	Recognizes chicken MHC Class II (B-L) molecules.
Reactivity:	Chicken
Purification:	Purified
Presentation:	PBS, sodium azide.
Usage Summary:	Flow cytometry. Immunoprecipitation. Western Blot. Immunohistochemistry. Flow Cytometry: 1 ug/1E6 cells.
Uses:	IHC, Western blot, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences