

MHC Class II I-A Rat anti-Mouse Monoclonal (FITC) Antibody - LS-C72635 - LSBio

CatalogID:	LS-C72635
Target:	MHC Class II I-A
Host	Rat
Clonality:	Monoclonal
Isotype:	IgG2b
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Mouse
Specificity:	Reacts specifically with a non-polymorphic I-A-encoded epitope on MHC Class II antigens.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, sodium azide.
Usage Summary:	Flow cytometry. Complement-dependent cytotoxicity assays. Flow Cytometry: 1 ug/1E6 cells.
Uses:	Flow Cytometry, Cytotoxic Assay (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences