

Monocytes, Granulocytes Mouse anti-Pig Monoclonal (FITC) Antibody - LS-C72757 - LSBio	
CatalogID:	LS-C72757
Target:	Monocytes, Granulocytes
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Pig
Specificity:	Reacts with >90% of porcine monocytes and granulocytes and < 5% of peripheral blood lymphocytes.
Reactivity:	Pig
Purification:	Purified
Presentation:	PBS, sodium azide.
Recommended Storage:	Short term: 4°C. Long term: Aliquot and add glycerol (40-50%). Store at -20°C. Avoid freeze-thaw cycles.
Usage Summary:	Flow cytometry. Immunoprecipitation. Flow Cytometry: 1 ug/1E6 cells.
Uses:	Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	