

Monocytes, Granulocytes Mouse anti-Pig Monoclonal (PE) Antibody - LS-C72759 - LSBio	
CatalogID:	LS-C72759
Target:	Monocytes, Granulocytes
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Pig
Specificity:	Reacts with >90% of porcine monocytes and granulocytes and < 5% of peripheral blood lymphocytes.
Reactivity:	Pig
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Usage Summary:	Flow cytometry. Immunoprecipitation. Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 10/2/2014 © 2014 LifeSpan BioSciences	