

Monocytes, Macrophages Mouse anti-Chicken Monoclonal (FITC) Antibody - LS-C72764 - LSBio	
CatalogID:	LS-C72764
Target:	Monocytes, Macrophages
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Chicken
Specificity:	Recognizes chicken monocytes and macrophages. Does not react with B (Bu-1+) or T (CD3+) lymphocytes.
Reactivity:	Chicken
Purification:	Purified
Presentation:	PBS, sodium azide.
Usage Summary:	Flow cytometry. Immunohistochemistry. Flow Cytometry: 1 ug/1E6 cells.
Uses:	IHC, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	