

**Phaseolus vulgaris L Lectin Rabbit anti-Phaseolus vulgaris Polyclonal (FITC) Antibody - LS-C74294
 - LSBio**

CatalogID:	LS-C74294
Target:	Phaseolus vulgaris L Lectin
Host	Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Phaseolus vulgaris
Antigen Type:	Purified protein
Immunogen:	Phaseolus vulgaris ul (PHA-L, Red Kidney Bean).
Specificity:	Recognizes Phaseolus vulgaris L.
Reactivity:	Phaseolus vulgaris
Purification:	Ammonium sulfate precipitation
Presentation:	PBS, pH 7.2, 0.1% sodium azide, 25% glycerol.
Usage Summary:	Suitable for use in Immunodiffusion and IEP.
Uses:	Immunodiffusion, Immunoelectrophoresis (Optimal dilution to be determined by the researcher)
Size:	1,000 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences