

Alkaline Phosphatase Mouse anti-Human Monoclonal (PE) Antibody - LS-C74831 - LSBio

CatalogID:	LS-C74831
Target:	Alkaline Phosphatase
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Human
Antigen Type:	Purified protein
Immunogen:	Human alkaline phosphatase.
Specificity:	Recognizes human Alkaline Phosphatase.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 0.5% BSA, 0.1% sodium azide.
Recommended Storage:	Maybe stored at 4°C before opening. DO NOT FREEZE! Stable at 4°C as an undiluted liquid. Stable for 12 months at 4°C. Sensitive to light
Usage Summary:	Suitable for use in Flow Cytometry.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Concentration:	0.03 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences