

Phosphotyrosine Mouse Monoclonal (FITC) Antibody - LS-C74927 - LSBio

CatalogID:	LS-C74927
Target:	Phosphotyrosine
Host	Phosphotyrosine antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b,k
Conjugations:	Fluorescein (FITC)
Antigen Type:	Hapten
Immunogen:	Phosphotyrosine antibody was raised against phosphotyramine-KLH.
Specificity:	Recognizes tyrosine-phosphorylated proteins from all species.
Purification:	Purified
Presentation:	PBS, 0.05% sodium azide.
Usage Summary:	Suitable for use in Immunocytochemistry. Immunocytochemistry: 1.5-15 ug/ml of this lot detected phosphotyrosine-containing proteins in EGF stimulated human A431 cells.
Uses:	ICC (1.5 - 15 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences