

MHC Class II RLA-DR Mouse anti-Rabbit Monoclonal (PE) Antibody - LS-C75494 - LSBio	
CatalogID:	LS-C75494
Target:	MHC Class II RLA-DR
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG2b,k
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Rabbit
Antigen Type:	Cells
Immunogen:	Rabbit leukemic B-cell line.
Specificity:	Recognizes rabbit RLA-DR chain (Mr ~35kD). Reacts with RLA-DR-transfected cells, but not with RLA-DQ-transfected cells and, in normal adult rabbits, 50-60% of peripheral blood leukocytes and splenocytes. It has been reported that the antibody stains: lymphocytes in the thymic medulla and some cells in the cortex; periportal cells in the liver; and synovial lining cell (i. e., infiltrating macrophages) in inflamed synovium.
Reactivity:	Rabbit
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Usage Summary:	Flow cytometry. Immunohistochemistry (acetone-fixed, frozen tissue sections). Immunoprecipitation. Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014 © 2014 LifeSpan BioSciences	