

**S-Tag Mouse Monoclonal (FITC) Antibody - LS-C75967 - LSBio**

<b>CatalogID:</b>	LS-C75967
<b>Target:</b>	S-Tag
<b>Host</b>	S-Tag antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a,k
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Antigen Type:</b>	Recombinant protein
<b>Immunogen:</b>	S-Tag antibody was raised against s-Tag recombinant protein.
<b>Specificity:</b>	Recognizes S-Tag.
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, sodium azide.
<b>Usage Summary:</b>	Flow cytometry. Enzyme-Linked-Immunosorbent-Assay (ELISA); Clone 5K159 can be paired with 5K. Protein purification. Flow Cytometry: 1 ug/1E6 cells.
<b>Uses:</b>	Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	500 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences