

S-Tag Goat Polyclonal (PE) Antibody - LS-C75977 - LSBio

CatalogID:	LS-C75977
Target:	S-Tag
Host	S-Tag antibody was produced in Goat
Clonality:	Polyclonal
Conjugations:	Phycoerythrin (PE)
Antigen Type:	Recombinant protein
Immunogen:	S-Tag antibody was raised against s-Tag recombinant protein.
Specificity:	Recognizes S-Tag.
Purification:	Affinity purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Usage Summary:	Flow cytometry; Detection for recombinant protein with S-Tag. Enzyme-Linked-Immunosorbent-Assay (ELISA); Detection for recombinant protein with S-Tag. Purification for recombinant protein with S-Tag. Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences