

Streptococcus A Goat anti-Streptococcus Polyclonal (FITC) Antibody - LS-C76085 - LSBio	
CatalogID:	LS-C76085
Target:	Streptococcus A
Host	Streptococcus A antibody was produced in Goat
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Streptococcus A antibody was raised against Streptococcus
Antigen Type:	Cells
Immunogen:	Streptococcus A antibody was raised against streptococci, group A.
Specificity:	Type specific carbohydrate of group A Streptococcus. Does not react with other Strep groups. Functions both as capture and reporter antibody. Conjugates well with colloidal gold.
Reactivity:	Streptococcus
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2, 0.1% sodium azide. No stabilizing proteins have been added.
Usage Summary:	Suitable for use in ELISA and Immunofluorescence. Immunofluorescence: 1:10-1:50. Acetone fixation of the antigen source is recommended prior to staining. ELISA: 1:20-1:200.
Uses:	Immunofluorescence (1:10 - 1:50), ELISA (1:20 - 1:200) (Optimal dilution to be determined by the researcher)
Size:	1,000 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	