

| TCR Alpha+Beta Hamster anti-Mouse Monoclonal (PE) (H57-597) Antibody - LS-C76266 - LSBio | |
|--|---|
| CatalogID: | LS-C76266 |
| Target: | T Cell Receptor Alpha/Beta |
| Host | T Cell Receptor Alpha/Beta antibody was produced in Hamster |
| Clonality: | Monoclonal |
| Isotype: | IgG |
| Clone Name: | H57-597 |
| Conjugations: | Phycoerythrin (PE) |
| Immunogen Species: | TCR Alpha+Beta antibody was raised against Mouse |
| Antigen Type: | Purified protein |
| Immunogen: | TCR Alpha+Beta antibody was raised against tCR affinity purified from mouse T-cell hybridoma DO-11.10. |
| Specificity: | Recognizes common epitope of the beta chain of the mouse alpha/beta T cell receptor (TCR) complex. Does not react with gamma, delta TCR-bearing cells. |
| Reactivity: | Mouse |
| Purification: | Purified |
| Presentation: | PBS, sodium azide, a stabilizing agent.5. |
| Usage Summary: | Flow cytometry. Immunohistochemistry (acetone-fixed, frozen sections and formalin-fixed, paraffin-embedded sections). Immunoprecipitation. In vitro stimulation of a/b TCR-expressing T cells. In vitro and in vivo depletion of a/b TCR-bearing cells. Flow Cytometry: 0.2 ug/1E6 cells. |
| Uses: | IHC - Paraffin, IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher) |
| Size: | 100 µg |
| Requested From: | Japan |
| Laboratory Reagent For In Vitro Research Use Only | |
| Not for resale without prior written consent from LifeSpan BioSciences, Inc. | |
| Created on 10/2/2014 | |
| © 2014 LifeSpan BioSciences | |