

TCR Gamma+TCR Delta Hamster anti-Mouse Monoclonal (FITC) (GL3) Antibody - LS-C76283 - LSBio	
CatalogID:	LS-C76283
Target:	T Cell Receptor Gamma/Delta
Host	Hamster
Clonality:	Monoclonal
Isotype:	IgG
Clone Name:	GL3
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Mouse
Antigen Type:	Cells
Immunogen:	C57BL/6J intra epithelial lymphocytes.
Specificity:	Reacts with mouse TCR gamma/delta. Species cross-reactivity: Not established.
Reactivity:	Mouse
Purification:	Protein A purified
Presentation:	0.1% sodium azide, 0.5% BSA.
Usage Summary:	Flow Cytometry: Neat.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	