

TCR Gamma+TCR Delta Mouse anti-Chicken Monoclonal (PE) Antibody - LS-C76290 - LSBio	
CatalogID:	LS-C76290
Target:	T Cell Receptor Gamma/Delta
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Chicken
Specificity:	Recognizes chicken TCR gamma, delta.
Reactivity:	Chicken
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Usage Summary:	Flow cytometry. Immunohistochemistry (acetone-fixed, frozen sections). Immunoprecipitation. In ovo depletion of gd T cells.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	