

Toxoplasma gondii Rabbit anti-Toxoplasma gondii Polyclonal (FITC) Antibody - LS-C76441 - LSBio	
CatalogID:	LS-C76441
Target:	Toxoplasma gondii
Host	Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Toxoplasma gondii
Antigen Type:	Cells
Immunogen:	RH strain of T. gondii.
Specificity:	Purified trophozoites.
Reactivity:	Toxoplasma gondii
Purification:	Purified
Presentation:	PBS, pH 7.2, BSA (10 mg/ml), 0.1% sodium azide.
Usage Summary:	Suitable for use in Immunofluorescence. Immunofluorescence: 1:10-1:50. Direct FA staining of target antigens in a permissive tissue culture system. Acetone fixation of the antigen source is recommended prior to staining.
Uses:	Immunofluorescence (1:10 - 1:50) (Optimal dilution to be determined by the researcher)
Size:	1,000 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 10/2/2014 © 2014 LifeSpan BioSciences	