

BrdU Sheep Polyclonal Antibody - LS-C77258 - LSBio

CatalogID:	LS-C77258
Target:	Bromodeoxyuridine (BrdU)
Host	Bromodeoxyuridine (BrdU) antibody was produced in Sheep
Clonality:	Polyclonal
Isotype:	IgG2b
Antigen Type:	Purified protein
Immunogen:	BrdU antibody was raised against highly purified Bromodeoxyuridine.
Specificity:	Recognizes bromodeoxyuridine at ~307kD. Very low cross-reactivity to 5-methyl cytosine.
Purification:	Protein G purified
Presentation:	PBS, pH 7.5, 40% glycerol. No preservative added.
Usage Summary:	Suitable for use in ELISA, Western Blot, Immunohistochemistry and Immunoprecipitation. Also available conjugated to horseradish peroxidase.
Uses:	IHC, Western blot, Immunoprecipitation, ELISA (1:250) (Optimal dilution to be determined by the researcher)
Size:	50 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences