

BrdU Mouse Monoclonal Antibody - LS-C77284 - LSBio

CatalogID:	LS-C77284
Target:	Bromodeoxyuridine (BrdU)
Host	Bromodeoxyuridine (BrdU) antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG
Antigen Type:	Hapten
Immunogen:	BrdU antibody was raised against bromodeoxyuridine.
Specificity:	Binds to bromodeoxyuridine and cross-reacts with iodouridine (10%). Does not cross-react with fluorodeoxyuridine, nor with any endogenous cellular components such as thymidine or uridine.
Purification:	Tissue culture supernatant
Presentation:	No preservative.
Usage Summary:	Immunocytochemistry Flow Cytometry.
Uses:	ICC, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences