

Campylobacter jejuni Rabbit anti-Campylobacter jejuni Polyclonal (FITC) Antibody - LS-C77637 - LSBio

CatalogID:	LS-C77637
Target:	Campylobacter jejuni
Host	Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Campylobacter jejuni
Antigen Type:	Cells
Immunogen:	ATCC strain 29428.
Specificity:	Whole cells. This antiserum has not been absorbed alpha may react with related microorganisms.
Reactivity:	Campylobacter jejuni
Purification:	Purified
Presentation:	PBS, pH 7.2, BSA (10 mg/ml), 0.1% sodium azide.
Usage Summary:	Suitable for use in Immunofluorescence. Immunofluorescence: 1:10-1:50. Direct FA staining of target antigens in a permissive tissue culture system. Acetone fixation of the antigen source is recommended prior to staining.
Uses:	Immunofluorescence (1:10 - 1:50) (Optimal dilution to be determined by the researcher)
Size:	1,000 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences