

CD3 Rat Monoclonal (FITC) Antibody - LS-C78155 - LSBio

CatalogID:	LS-C78155
Target:	CD3
Host	CD3 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a
Conjugations:	Fluorescein (FITC)
Specificity:	Recognizes the mouse CD3 antigen which is expressed by mature T cells. NB. For optimal staining incubations should be performed at room temperature.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, pH 7.2, 0.1% sodium azide.
Usage Summary:	Flow Cytometry: Neat-1:10.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences