

**CD3 Mouse Monoclonal (PE) Antibody - LS-C78198 - LSBio**

<b>CatalogID:</b>	LS-C78198
<b>Target:</b>	CD3
<b>Host</b>	CD3 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Specificity:</b>	Recognizes CD3, a 19 kD cell surface protein which is expressed on >95% of circulating human peripheral blood T-cells. Species cross Reactivity: chimpanzee.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein G purified
<b>Reconstitution:</b>	deionized water
<b>Presentation:</b>	Lyophilized from PBS, pH 7.4, 1% BSA, 0.09% sodium azide. Phycoerythrin (RPE).
<b>Usage Summary:</b>	Suitable for use in Flow Cytometry. Flow Cytometry: Neat, Use 10 ul to label 1E6 cells in 100ul.
<b>Uses:</b>	Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 tst
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences