

SSEA-1 / Lewis x / CD15 Mouse Monoclonal (FITC) Antibody - LS-C78322 - LSBio

CatalogID:	LS-C78322
Target:	Stage-Specific Embryonic Antigen-1/Lewis X/CD15
Host	Mouse
Clonality:	Monoclonal
Isotype:	IgM
Conjugations:	Fluorescein (FITC)
Specificity:	Recognizes human CD15.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, sodium azide and a stabilizing agent
Usage Summary:	Flow cytometry. Immunohistochemistry (frozen sections). Immunoprecipitation. Inhibition of neutrophil phagocytosis and bactericidal activity. Flow Cytometry: 10 ul/1E6 cells.
Uses:	IHC - Frozen, Immunoprecipitation, Flow Cytometry, Functional Assay (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences