

CD16 + CD32 Rat Monoclonal (PE) (93) Antibody - LS-C78366 - LSBio	
CatalogID:	LS-C78366
Target:	CD16/CD32 (FcγR3/FcγR2)
Host	CD16/CD32 (FcγR3/FcγR2) antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2b,k
Clone Name:	93
Conjugations:	Phycoerythrin (PE)
Specificity:	Recognizes murine CD16/32 (CD16/Fc-gammall and CD32/Fc-gIII receptors). Recognizes a conformational epitope formed by Fc-gII and Fc-gIII receptors.
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Usage Summary:	Flow cytometry. Blocking of Fc-g receptors. Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	