

CD16 + CD32 Rat Monoclonal (PE) (93) Antibody - LS-C78370 - LSBio	
<b>CatalogID:</b>	LS-C78370
<b>Target:</b>	CD16/CD32 (FcγR3/FcγR2)
<b>Host</b>	CD16/CD32 (FcγR3/FcγR2) antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2b,k
<b>Clone Name:</b>	93
<b>Conjugations:</b>	Phycoerythrin (PE)
<b>Specificity:</b>	Recognizes murine CD16/32 (CD16/Fc-gamma1 and CD32/Fc-gamma3 receptors). Recognizes a conformational epitope formed by Fc-gamma1 and Fc-gamma3 receptors.
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Purified
<b>Presentation:</b>	PBS, 0.1% sodium azide, stabilizing agent.
<b>Usage Summary:</b>	Flow cytometry. Blocking of Fc-g receptors. Flow Cytometry: 0.2 ug/1E6 cells.
<b>Uses:</b>	Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	