

IL4 Rat anti-Mouse Monoclonal (PE) Antibody - LS-C7872 - LSBio

CatalogID:	LS-C7872
Target:	interleukin 4 (IL4)
Synonyms:	IL4 Antibody, B cell growth factor 1 Antibody, B-cell stimulatory factor 1 Antibody, B_cell stimulatory factor 1 Antibody, BCGF1 Antibody, BSF-1 Antibody, Binetrakin Antibody, BCGF-1 Antibody, Interleukin 4 Antibody, Interleukin-4 Antibody, Pitrakinra Antibody, BSF1 Antibody, IL-4 Antibody
Family / Subfamily:	Hormone / not assigned-Hormone
Host	IL4 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	IL4 antibody was raised against Mouse
Antigen Type:	Recombinant protein
Immunogen:	IL4 antibody was raised against recombinant murine IL-4.
Specificity:	Recognizes mouse IL-4.
Reactivity:	Mouse
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.1% sodium azide.
Recommended Storage:	Store at 4°C, do not freeze. Freezing R-Phycoerythrin conjugates will result in substantial loss of activity.
Usage Summary:	Suitable for use in Flow Cytometry. Flow Cytometry: Neat 10 ul labels 1E6 cells in 100ul. Membrane permeabilization is required for this application.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	500 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences