

Collagen IV Mouse Monoclonal (PE) Antibody - LS-C79598 - LSBio

CatalogID:	LS-C79598
Target:	Collagen IV
Host	Collagen IV antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE)
Specificity:	Reacts with human Type IV Collagen.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Store at +4°C. Do not freeze. Product is photosensitive and should be protected from light.
Usage Summary:	Immunohistochemistry (frozen sections). Enzyme-Linked-Immunosorbent-Assay (ELISA).
Uses:	IHC - Frozen, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Concentration:	0.5 mg/ml
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences