

Coronavirus OC-43 Nucleoprotein Mouse anti-Coronavirus Monoclonal Antibody - LS-C79764 - LSBio

<b>CatalogID:</b>	LS-C79764
<b>Target:</b>	Coronavirus OC-43 Nucleoprotein
<b>Host</b>	Mouse
<b>Clonality:</b>	Monoclonal
<b>Immunogen Species:</b>	Coronavirus
<b>Antigen Type:</b>	Purified protein
<b>Immunogen:</b>	Coronavirus Strain OC-43 nucleoprotein.
<b>Specificity:</b>	Recognizes Coronavirus Strain OC-43 nucleoprotein.
<b>Reactivity:</b>	Coronavirus
<b>Purification:</b>	Protein A/G purified
<b>Presentation:</b>	PBS, pH 7.6, 0.1% sodium azide
<b>Usage Summary:</b>	Suitable for use in ELISA.
<b>Uses:</b>	ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences