

IL13 Mouse anti-Human Monoclonal (FITC) Antibody - LS-C8140 - LSBio	
CatalogID:	LS-C8140
Target:	interleukin 13 (IL13)
Synonyms:	IL13 Antibody, Allergic rhinitis Antibody, ALRH Antibody, BHR1 Antibody, Interleukin-13 Antibody, NC30 Antibody, p600 Antibody, IL-13 Antibody, Interleukin 13 Antibody
Family / Subfamily:	Hormone / not assigned-Hormone
Host	IL13 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Fluorescein (FITC)
Immunogen Species:	IL13 antibody was raised against Human
Specificity:	Recognizes human IL-13. Blocks IL-13 induced proliferation.
Reactivity:	Human
Purification:	Purified
Presentation:	50 mM Tris, 100 mM sodium chloride, 0.02% sodium azide, pH 7.4, 1% BSA.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Suitable for use in Functional Studies. Functional Studies: Blocks IL-13 induced proliferation.
Uses:	Functional Assay (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	