

BrdU Mouse Monoclonal Antibody - LS-C83001 - LSBio

CatalogID:	LS-C83001
Target:	Bromodeoxyuridine (BrdU)
Host	Bromodeoxyuridine (BrdU) antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Antigen Type:	Hapten
Immunogen:	BrdU antibody was raised against bromodeoxyuridine-BSA conjugate.
Specificity:	Specifically reacts with bromodeoxyuridine.
Presentation:	Lyophilized from phosphate buffered saline, pH 7.4, 0.1% sodium azide
Recommended Storage:	Short term 4°C, long term aliquot and store at -20°C, avoid freeze thaw cycles.
Usage Summary:	Immunoassay: 5 ug/ml.
Uses:	Undetermined
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences