

Myc Tag Rabbit anti-Human Polyclonal Antibody - LS-C83258 - LSBio

CatalogID:	LS-C83258
Target:	Myc Tag
Host	Myc Tag antibody was produced in Rabbit
Clonality:	Polyclonal
Immunogen Species:	Myc Tag antibody was raised against Human
Immunogen:	Myc Tag antibody was raised against c-myc (EQKLISEEDL) coupled to KLH. Percent identity by BLAST analysis: Human, Chimpanzee, Gorilla, Mouse, Rat (100%); Orangutan, Monkey, Marmoset, Tamarin, Panda, Cat (90%); Horse, Pufferfish, Paramecium, Salpingoeca (80%).
Reactivity:	Human, Chimpanzee, Gorilla, Mouse, Rat
Predicted Reactivity:	Orangutan, Monkey
Presentation:	PBS, pH 7.2, 0.1% sodium azide
Recommended Storage:	Store at 4°C, stable for one year.
Usage Summary:	ELISA: Coating: 2-10 ug/ml, Primary: 0.03-1 ug/m. Western Blot: Colorimetric: 0.1-1 ug/ml Chemiluminescent: 0.03-1 ug/ml.
Uses:	IHC, ICC, Western blot, ELISA (0.03 - 2 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences