

GST / Glutathione S-Transferase Goat anti-Schistosoma japonicum Polyclonal (FITC) Antibody - LS-C83273 - LSBio	
CatalogID:	LS-C83273
Target:	GST / Glutathione S-Transferase
Host	Goat
Clonality:	Polyclonal
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Schistosoma japonicum
Antigen Type:	Purified protein
Immunogen:	GST from Schistosoma japonicum coupled to KLH.
Reactivity:	Schistosoma japonicum
Presentation:	PBS, pH 7.2, 0.1% sodium azide
Recommended Storage:	Store at 4°C, stable for one year.
Uses:	IHC, ICC (0.5 - 5 µg/ml) (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	