

KLRG1 Hamster anti-Mouse Monoclonal (PE) Antibody - LS-C8380 - LSBio	
CatalogID:	LS-C8380
Target:	killer cell lectin-like receptor subfamily G, member 1 (KLRG1)
Synonyms:	KLRG1 Antibody, 2F1 Antibody, CLEC15A Antibody, MAFA-2F1 Antibody, MAFA-L Antibody, MAFA-LIKE Antibody, MAFA-like receptor Antibody, MAFAL Antibody
Host	KLRG1 antibody was produced in Hamster
Clonality:	Monoclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	KLRG1 antibody was raised against Mouse
Antigen Type:	Cells
Immunogen:	KLRG1 antibody was raised against activated NK (A-LAK) cells from B6 (H-2b) mice.
Specificity:	Reacts with mouse killer cell lectin-like receptor G1 (KLRG1, Mr 30-38kD).
Reactivity:	Mouse
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	Long term: Add glycerol (40-50%) -20°C; Short term: +4°C
Usage Summary:	Suitable for use in Flow Cytometry and Immunoprecipitation. Flow Cytometry: 0.2 ug/1E6 cells.
Uses:	Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	