

**BrdU Mouse Monoclonal (BU5.1) Antibody - LS-C84809 - LSBio**

<b>CatalogID:</b>	LS-C84809
<b>Target:</b>	Bromodeoxyuridine (BrdU)
<b>Host</b>	Bromodeoxyuridine (BrdU) antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Clone Name:</b>	BU5.1
<b>Antigen Type:</b>	Hapten
<b>Immunogen:</b>	BrdU antibody was raised against bromodeoxyuridine.
<b>Specificity:</b>	This antibody is specific for bromodeoxyuridine and does not cross react with thymidine.
<b>Purification:</b>	Protein A purified
<b>Reconstitution:</b>	PBS
<b>Presentation:</b>	Lyophilized from PBS, pH 7.4, 0.5% BSA, 0.09 % sodium azide
<b>Recommended Storage:</b>	Prior to reconstitution, +4°C. Following reconstitution store at -20°C. Avoid freeze-thaw cycles.
<b>Uses:</b>	Undetermined
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences