

THY1 / CD90 Mouse anti-Human Monoclonal (F15-42-1) Antibody - LS-C85050 - LSBio

CatalogID:	LS-C85050
Target:	Thy-1 cell surface antigen (THY1)
Synonyms:	THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody
Host	THY1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	F15-42-1
Immunogen Species:	THY1 / CD90 antibody was raised against Human
Specificity:	Thy-1 molecule is heavily glycosylated membrane bound protein of 18kD.
Reactivity:	Human
Purification:	Ion exchange chromatography
Presentation:	PBS, 10 mM sodium azide, 1 mg/ml BSA
Recommended Storage:	Store at 4°C for short term applications. For long term storage, aliquot and store at -20°C.
Uses:	Undetermined
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences