

THY1 / CD90 Rat anti-Mouse Monoclonal (G7) Antibody - LS-C85051 - LSBio

CatalogID:	LS-C85051
Target:	Thy-1 cell surface antigen (THY1)
Synonyms:	THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody
Host	THY1 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2c,k
Clone Name:	G7
Immunogen Species:	THY1 / CD90 antibody was raised against Mouse
Immunogen:	THY1 / CD90 antibody was raised against mouse T-Cell hybridoma c6/G8, produced by fusion of porcine insulin-specific BALB/c T cells with the AKR thymoma line BW 5147.
Specificity:	This antibody reacts with the Thy-1.1 and Thy- 1.2 alloantigens on thymocytes, peripheral T Lymphocytes, and some interepithelial T lymphocytes (IEL, DEC) and neurons, but not B lymphocytes, of all mouse strains.
Reactivity:	Mouse
Presentation:	PBS
Recommended Storage:	Store at 4°C. Do not freeze.
Uses:	Undetermined
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences