

THY1 / CD90 Rat anti-Mouse Monoclonal (30-H12) Antibody - LS-C85052 - LSBio

CatalogID:	LS-C85052
Target:	Thy-1 cell surface antigen (THY1)
Synonyms:	THY1 Antibody, CD90 antigen Antibody, CD90 Antibody, CDw90 Antibody, Thy-1 antigen Antibody, Thy-1 T-cell antigen Antibody, Thy-1 Antibody, Thy-1 cell surface antigen Antibody, Thy-1 membrane glycoprotein Antibody
Host	THY1 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2b,k
Clone Name:	30-H12
Immunogen Species:	THY1 / CD90 antibody was raised against Mouse
Specificity:	This antibody does not cross react with mouse and rat strains bearing the Thy-1.1 alloantigen.
Reactivity:	Mouse
Presentation:	100 mM BBS, pH 8.2
Recommended Storage:	Store at 4°C. Do not freeze.
Uses:	IHC, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences