

**BrdU Sheep Polyclonal Antibody - LS-C85687 - LSBio**

<b>CatalogID:</b>	LS-C85687
<b>Target:</b>	Bromodeoxyuridine (BrdU)
<b>Host</b>	Bromodeoxyuridine (BrdU) antibody was produced in Sheep
<b>Clonality:</b>	Polyclonal
<b>Immunogen:</b>	BrdU antibody was raised against highly pure Bromodeoxyuridine.
<b>Specificity:</b>	Not species specific. Nucleotide is present in all species. Very low cross reactivity to 5-methyl cytosine.
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, no preservative
<b>Recommended Storage:</b>	Long term: -20°C; Short term: +4°C; Avoid freeze-thaw cycles.
<b>Uses:</b>	IHC, Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	250 µg
<b>Concentration:</b>	1.93 mg/ml
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences