

LAMP1 / CD107a Mouse anti-Human Monoclonal (PE) Antibody - LS-C8577 - LSBio

CatalogID:	LS-C8577
Target:	lysosomal-associated membrane protein 1 (LAMP1)
Synonyms:	LAMP1 Antibody, CD107a Antibody, CD107a antigen Antibody, LAMP-1 Antibody, LGP120 Antibody, LAMPA Antibody
Host	LAMP1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	LAMP1 / CD107a antibody was raised against Human
Specificity:	Recognizes human LAMP 1.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	+4°C, avoid freezing
Usage Summary:	Suitable for use in Flow Cytometry and Immunohistochemistry. Flow Cytometry: 10 ul/1E6 cells. Immunohistochemistry: acetone-fixed, frozen tissue sections.
Uses:	IHC - Frozen, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences