

FITC Rabbit Polyclonal Antibody - LS-C85811 - LSBio

CatalogID:	LS-C85811
Target:	Fluorescein Isothiocyanate (FITC)
Host:	Rabbit
Clonality:	Polyclonal
Antigen Species:	
Antigen Type:	Hapten
Immunogen:	Fluorescein isothiocyanate-KLH.
Specificity:	Fluorescein 100%
Purification:	Antiserum
Presentation:	Neat serum, no preservatives added
Recommended Storage:	Long term: -20°C; Short term: +4°C; Avoid freeze-thaw cycles.
Uses:	Immunoprecipitation (IP) (Optimal dilution to be determined by the researcher)
Size:	1,000 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 3/16/2012

© 2012 LifeSpan BioSciences