

LAMP2 / CD107b Mouse anti-Human Monoclonal (PE) Antibody - LS-C8594 - LSBio

CatalogID:	LS-C8594
Target:	lysosomal-associated membrane protein 2 (LAMP2)
Synonyms:	LAMP2 Antibody, CD107b Antibody, CD107b antigen Antibody, LAMPB Antibody, LGP110 Antibody, LAMP-2 Antibody
Family / Subfamily:	Cancer
Host	LAMP2 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	LAMP2 / CD107b antibody was raised against Human
Specificity:	Recognizes human LAMP 2.
Reactivity:	Human
Purification:	Purified
Presentation:	PBS, 0.1% sodium azide, stabilizing agent.
Recommended Storage:	+4°C, avoid freezing
Usage Summary:	Suitable for use in Flow Cytometry and Immunohistochemistry. Flow Cytometry: 10 ul/1E6 cells. Immunohistochemistry: Acetone-fixed, frozen tissue sections.
Uses:	IHC - Frozen, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences