

IAPP / Amylin Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C85977 - LSBio

CatalogID:	LS-C85977
Target:	islet amyloid polypeptide (IAPP)
Synonyms:	IAPP Antibody, Diabetes-associated peptide Antibody, Islet amyloid polypeptide Antibody, AMYLIN Antibody, IAP Antibody, Insulinoma amyloid peptide Antibody
Host	IAPP antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	Fluorescein (FITC)
Immunogen Species:	IAPP / Amylin antibody was raised against Human
Immunogen:	IAPP / Amylin antibody was raised against amylin (# 7321) from Peninsula.
Specificity:	Human, rat, cat, pig and guinea pig.
Reactivity:	Human
Reconstitution:	PBS. Possible additional volumes for resuspension: 50 µl
Presentation:	Lyophilized antiserum
Recommended Storage:	Prior to reconstitution, +4°C. Following reconstitution store at -20°C. Avoid freeze-thaw cycles.
Uses:	IHC (Optimal dilution to be determined by the researcher)
Size:	50 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences