

Rat IgG Donkey anti-Rat Polyclonal (Rhodamine) Antibody - LS-C86444 - LSBio

CatalogID:	LS-C86444
Target:	Rat IgG
Host	Rat IgG antibody was produced in Donkey
Clonality:	Polyclonal
Conjugations:	Rhodamine
Immunogen Species:	Rat IgG antibody was raised against Rat
Specificity:	Minimal Crossreactivity with Bovine, Chicken, Goat, Guinea Pig, Syrian Hamster, horse, human, rabbit and sheep
Reactivity:	Rat
Purification:	Affinity purified
Reconstitution:	Distilled water
Presentation:	Lyophilized from 0.01 M sodium phosphate, 0.25 M sodium chloride, pH 7.6, 15 mg/ml BSA, 0.05% sodium azide
Recommended Storage:	Store vial at 4 C prior to opening. This product is stable at 4 C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, mix with an equal volume of glycerol, aliquot contents and freeze at -20 C or below.
Uses:	Undetermined
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences