

Chicken IgY Donkey anti-Chicken Polyclonal (Rhodamine) Antibody - LS-C86449 - LSBio

CatalogID:	LS-C86449
Target:	Chicken IgY
Host	Chicken IgY antibody was produced in Donkey
Clonality:	Polyclonal
Conjugations:	Rhodamine
Immunogen Species:	Chicken IgY antibody was raised against Chicken
Specificity:	Minimal Cross to Bovine,Goat,Guinea Pig,Syrian Hamster, Horse, Human, Mouse, Rat, rabbit and Sheep Serum Proteins.
Reactivity:	Chicken
Purification:	Affinity purified
Reconstitution:	Distilled water
Presentation:	Lyophilized from 0.01 M sodium phosphate, 0.25 M sodium chloride, pH 7.6, 15 mg/ml BSA, 0.05% sodium azide
Recommended Storage:	Prior to reconstitution, +4°C. Following reconstitution store at -20°C. Avoid freeze-thaw cycles.
Uses:	Undetermined
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences