

E. coli O157:H7 Goat anti-E. coli Polyclonal (FITC) Antibody - LS-C86585 - LSBio

CatalogID:	LS-C86585
Target:	E. coli O157:H7
Host	E. coli O157:H7 antibody was produced in Goat
Clonality:	Polyclonal
Conjugations:	Fluorescein (FITC)
Immunogen Species:	E. coli O157:H7 antibody was raised against E. coli
Antigen Type:	Cells
Immunogen:	E. coli O157:H7 antibody was raised against e. coli serotype O157: H7.
Specificity:	This antibody is highly specific for E coli 0157:H7.
Reactivity:	E. coli
Purification:	Affinity purified
Reconstitution:	Water
Presentation:	Lyophilized
Recommended Storage:	4°C while lyophilized, once reconstituted aliquot and freeze at -20°C
Uses:	Undetermined
Size:	500 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences