

MYC / c-Myc Goat anti-Human Polyclonal (aa410-419) (HRP) Antibody - LS-C86627 - LSBio

CatalogID:	LS-C86627
Target:	v-myc avian myelocytomatosis viral oncogene homolog
Synonyms:	MYC Antibody, BHLHe39 Antibody, C-Myc Antibody, MRTL Antibody, Myc proto-oncogene protein Antibody, Proto-oncogene c-Myc Antibody, Rats1 Antibody, Transcription factor p64 Antibody, v-myc Antibody
Host	MYC antibody was produced in Goat
Clonality:	Polyclonal
Conjugations:	Horseradish Peroxidase (HRP)
Immunogen Species:	MYC / c-Myc antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MYC / c-Myc antibody was raised against synthetic peptide rep. amino acid residues 410-419 of human myc conjugated to.
Epitope:	aa410-419
Reactivity:	Human
Purification:	Affinity purified
Presentation:	50 mM HEPES, pH 7.1, 0.1 M sodium chloride, 0.2% BSA, 0.01% Thimerosal
Recommended Storage:	Store at 4°C. Do not freeze.
Uses:	Undetermined
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences