

Streptococcus A Rabbit anti-Streptococcus Polyclonal (FITC) Antibody - LS-C86701 - LSBio

CatalogID:	LS-C86701
Target:	Streptococcus A
Host	Streptococcus A antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Streptococcus A antibody was raised against Streptococcus
Immunogen:	Streptococcus A antibody was raised against streptococci Group A (STREP A).
Specificity:	Streptococci Group A (STREP A) 100%
Reactivity:	Streptococcus
Presentation:	PBS, pH 7.2, 1% BSA, 0.1% sodium azide
Recommended Storage:	Long term: -20°C; Short term: +4°C; Avoid freeze-thaw cycles.
Uses:	Undetermined
Size:	1,000 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/2/2014

© 2014 LifeSpan BioSciences