

Streptococcus B Rabbit anti-Streptococcus Polyclonal (FITC) Antibody - LS-C86703 - LSBio	
CatalogID:	LS-C86703
Target:	Streptococcus B
Host	Streptococcus B antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	Fluorescein (FITC)
Immunogen Species:	Streptococcus B antibody was raised against Streptococcus
Immunogen:	Streptococcus B antibody was raised against streptococci Group B (STREP B).
Specificity:	Streptococci Group B (Strep. B) 100%
Reactivity:	Streptococcus
Purification:	Column chromatography
Presentation:	PBS, pH 7.2, 1% BSA, 0.1% sodium azide
Recommended Storage:	Long term: -20°C; Short term: +4°C; Avoid freeze-thaw cycles.
Uses:	Immunofluorescence (1:20 - 1:100) (Optimal dilution to be determined by the researcher)
Size:	1,000 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	