

LAMP1 / CD107a Mouse anti-Human Monoclonal (FITC) (H4A3) Antibody - LS-C86764 - LSBio	
CatalogID:	LS-C86764
Target:	lysosomal-associated membrane protein 1 (LAMP1)
Synonyms:	LAMP1 Antibody, CD107a Antibody, CD107a antigen Antibody, LAMP-1 Antibody, LGP120 Antibody, LAMPA Antibody
Host	LAMP1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	H4A3
Conjugations:	Fluorescein (FITC)
Immunogen Species:	LAMP1 / CD107a antibody was raised against Human
Specificity:	Specific for Human CD107a/LAMP-1
Reactivity:	Human
Presentation:	PBS, sodium azide
Recommended Storage:	Store at 4°C. Do not freeze.
Uses:	Undetermined
Size:	100 tst
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 10/2/2014	
© 2014 LifeSpan BioSciences	